



CADD Pump Instructions



1. Clean and disinfect your workspace.
2. Take out all your supplies: medication, tubing, flushes, and alcohol wipes.
3. Wash your hands.
4. Stop the pump and disconnect the tubing from your IV line.
 - a. The batteries should be changed every 3 days. Now is the best time to change them. Turn off the pump, change the batteries (found on top of the pump), and turn the pump back on.
5. Remove the port protector on the IV bag.
6. Open the tubing, remove the white cap protecting the spike, and attach the administration set to the bag.
7. Remove the blue tab on the administration set "cassette."
8. Remove the old tubing from the pump (cassette) and attach the new one. (when looking at the pump, on the right-hand side, pull the silver lever towards you).
9. Attach the new tubing (cassette) by inserting the hooks on the cassette into the hinge pins on the bottom left of the pump. Swing the cassette up and then push the silver lever closed.
10. After you attach the new tubing, you will receive two messages on your screen (reset reservoir volume and prime tubing). The answer should be "YES" to both questions.
11. Swab the end of your IV line with alcohol for 20 to 30 seconds and let dry.
12. Flush your line with a normal saline syringe flush.
13. Connect the CADD tubing to your IV line.
14. Press "Start" on the pump.
15. Place medication and CADD pump into the carrying bag.
16. Clean your workspace.
17. Wash your hands.



- This information is for reference only and should only be performed after being taught by your nurse.
- Please call us if you have any questions or experience: redness, swelling at the IV site, unable to flush or infuse your medication, IV becomes dislodged, IV dressing has fallen off, you have a temperature over 101, you begin to experience any negative side effects of the medication.
- *****Call 911 if you are having shortness of breath or chest pain*****