




### Mini-Bag Plus Instructions

1. Clean and disinfect your workspace.
2. Take out all your supplies: medication, tubing, flushes, and alcohol wipes.
3. Mix the medication.
  - a. Bend up, then down to break the seal.
  - b. Hold the bag with the vial down. Squeeze the solution into the vial until half full. Shake gently to mix the powder.
  - c. Hold the bag with the vial upside down. Squeeze the bag to force air into the vial. Release the bag, and the medication should drain out of the vial and into the bag.
  - d. Ensure the medication is completely dissolved. Repeat steps B and C until the vial is empty of the medication and the solution is thoroughly mixed. Leave the vial attached.



4. Wash your hands.
  5. Swab the end of your IV line with alcohol for 20 to 30 seconds and let dry
  6. Flush your line with a normal saline syringe flush
  7. Remove the port protector and attach the administration set.
  8. Hang your bag on the IV pole, unclamp the clamp on the administration set, prime the admin set, set the flow rate as indicated on the label, connect the admin set to your injection site, and infuse the medication.
  9. When the bag is empty, disconnect the admin set from yourself.
  10. Flush your line as directed by your nurse.
  11. Clean your workspace.
  12. Wash your hands
- This information is for reference only and should only be performed after being taught by your nurse.
  - Please call us if you have any questions or experience: redness, swelling at the IV site, unable to flush or infuse your medication, IV becomes dislodged, IV dressing has fallen off, you have a temperature over 101, you begin to experience any negative side effects of the medication.
  - **\*\*\*Call 911 if you are having shortness of breath or chest pain\*\*\***